

Air/Oxygen Blenders

Blender Gallery

Blender Gallery

Optimize Patient Care without Compromising on Quality



Part Number: 2003FF

Flow Range: 0-30 lpm

Output Port Flow Range:

Left: 0-15 lpm Flowmeter (Green O₂ Label)

Right: 0-3.5 lpm Flowmeter (Green O₂ Label)

Bleed flow activated and deactivated by rotating the flowmeter

Bleed for Flows <3 lpm: 3 lpm

FiO₂ Accuracy with Bleed: 0-3 lpm

FiO₂ Accuracy without Bleed: 3-30 lpm

**2003FF
LOW FLOW
BLENDER**

Oxygen Administration - Bubble CPAP Therapy - HFNC Therapy

Our Bio-Med Devices Blenders offer:

- A full range of high-quality, precision-engineered configurations
- Accurate oxygen therapy while maximizing cost-efficiency
- Top-tier craftsmanship that ensures reliable performance for your most vulnerable patients
- Industry-leading durability and unmatched support from our dedicated Service Team

MANUFACTURED BY:



Bio-Med Devices, Inc.
61 Soundview Road
Guilford, CT 06437

Proudly Made in the USA

Telephone: 1-800-224-6633 or 203-458-0202

Fax: 203-458-0440


Website: www.biomeddevices.com

Email: custserv@biomeddevices.com

DDOC063 Rev031925

Air/Oxygen Blenders

Blender Gallery



**Low Flow
2003**

Part Number: 2003

Flow Range: 0-30 lpm

Output Port Flow Range:

Left: 0-30 lpm

Right: Bleed flow activated when flowmeter is engaged

Bleed for Flows <3 lpm: 3 lpm



**Low Flow
2003M**

Part Number: 2003M (MR Conditional)

Flow Range: 0-30 lpm

Output Port Flow Range:

Left: 0-30 lpm

Right: Bleed flow activated when flowmeter is engaged

Bleed for Flows <3 lpm: 3 lpm



**Low Flow
2003F**

Part Number: 2003F

Flow Range: 0-30 lpm

Output Port Flow Range:

Left: Plugged

Right: 0-15 lpm Flowmeter (Green O₂ Label)

Bleed flow activated and deactivated by rotating the flowmeter

Bleed for Flows <3 lpm: 3 lpm

FiO₂ Accuracy with Bleed: 0-30 lpm

FiO₂ Accuracy without Bleed: 3-30 lpm



**Low Flow
2003FL**

Part Number: 2003FL

Flow Range: 0-30 lpm

Output Port Flow Range:

Left: 0-15 lpm Flowmeter (Green O₂ Label)

Bottom: 0-30 lpm

Right: Knob to switch bleed on & off

Bleed for Flows <3 lpm: 3 lpm

FiO₂ Accuracy with Bleed: 0-30 lpm

FiO₂ Accuracy without Bleed: 3-30 lpm

MANUFACTURED BY:



Bio-Med Devices, Inc.
61 Soundview Road
Guilford, CT 06437

Proudly Made in the USA

Telephone: 1-800-224-6633 or 203-458-0202

Fax: 203-458-0440

Website: www.biomeddevices.com

Email: custserv@biomeddevices.com



BIO-MED
DEVICES

Air/Oxygen Blenders

Blender Gallery



Low Flow
2003P

Part Number: 2003P

Flow Range: 0-30 lpm

Output Port Flow Range:

Left: 0-30 lpm

Bottom: 0-30 lpm

Right: Bleed flow activated when flowmeter is engaged

Bleed for Flows <3 LPM: 3 lpm



High-Low Flow
2004

Part Number: 2004

Flow Range: 0-120 lpm

Output Port Flow Range:

Left: 0-120 lpm

Bottom: 0-120 lpm

Right: Bleed flow activated when flowmeter is engaged

Bleed for Flows <15 lpm: 12 lpm



High-Low Flow
2004F7015

Part Number: 2004F7015

Flow Range: 0-120 lpm

Output Port Flow Range:

Left: 0-70 lpm Flowmeter (Green O₂ Label)

Bottom: 0-120 lpm

Right: 0-15 lpm Flowmeter (Green O₂ Label)

Bleed for Flows <15 lpm: 12 lpm

FiO₂ Accuracy with Bleed: 2-108 lpm

FiO₂ Accuracy without Bleed: 15-120 lpm



High-Low Flow
2001

Part Number: 2001

Flow Range: 0-120 lpm

Output Port Flow Range:

Left: Plugged

Bottom: 0-120 lpm

Right: Bleed flow activated when flowmeter is engaged

Bleed for Flows <15 lpm: 12 lpm

MANUFACTURED BY:

Bio-Med Devices, Inc.
61 Soundview Road
Guilford, CT 06437

Proudly Made in the USA

Telephone: 1-800-224-6633 or 203-458-0202

Fax: 203-458-0440

Website: www.biomeddevices.com

Email: custserv@biomeddevices.com

DDOC063 Rev031925

Air/Oxygen Blenders

Blender Gallery



High-Low Flow 2001K

Part Number: 2001K

Flow Range: 0-120 lpm

Output Port Flow Range:

Left: Plugged

Bottom: 0-120 lpm

Right: Knob to switch bleed on & off

Bleed for Flows <15 lpm: 12 lpm

FiO2 Accuracy with Bleed: 2-108 lpm

FiO2 Accuracy without Bleed: 15-120 lpm



High Flow 2002

Part Number: 2002

Flow Range: 0-120 lpm

Output Port Flow Range:


Left: Plugged

Bottom: 0-120 lpm

Right: Plugged

Bleed: None

FiO2 Accuracy: 15-120 lpm



High Flow 2002F70D

Part Number: 2002F70D

Flow Range: 0-120 lpm

Output Port Flow Range:

Left: 0-70 lpm Flowmeter (Green O₂ Label)

Bottom: 0-120 lpm

Right: Plugged

Bleed: None

FiO2 Accuracy: 15-120 lpm



High Flow 2002F70

Part Number: 2002F70

Flow Range: 0-120 lpm

Output Port Flow Range:

Left: 0-70 lpm Flowmeter (White O₂ Label)

Right: Plugged

Bleed: None

FiO2 Accuracy: 15-120 lpm

MANUFACTURED BY:



Bio-Med Devices, Inc.
61 Soundview Road
Guilford, CT 06437

Proudly Made in the USA

Telephone: 1-800-224-6633 or 203-458-0202

Fax: 203-458-0440

Website: www.biomeddevices.com

Email: custserv@biomeddevices.com

DDOC063 Rev031925