

TALENT II

DTP886

Specifications



Configurations

No.	Brand	Model	Product Name	Quantity
1	Angell	Angell-DR-HCP	Dynamic Flat Panel Detector	1
2	Angell	TS-5	DXRay Diagnostic Workstation	1
3	Angell	XDRP80	High Frequency Generator	1
4	Angell/ Siemens	SV150/40	X-Ray Tube	1
5	Angell	XS-F	Collimator	1
6	Angell	DC-DT-4	Remote Control Table	1
7	Angell	YKT-4	Remote-Control Console	1
8	Angell	YKT-4M	Movable Console with Monitor (Optional)	1
9	JPI	JPI GRID-1000	Grid	1
10	/	/	Report Printer	1
11	/	/	Microphone/Speaker	1
12	/	/	High Voltage Cables	2
13	/	/	TOMO Wide range linear scanning depth Imaging (Optional)	1

High Frequency Generator

Specifications

- Model: XDRP80
- Brand: Angell
- Output power: 80kW
- Input Power Frequency: 50Hz/60Hz
- Output Voltage for Radiography: 40kV-150kV
- Output Voltage for Fluoroscopy: 40kV-125kV
- Output Current for Radiography: 10mA-1000mA
- Output Current for Fluoroscopy: 0.5mA-20mA
- mAs Range: 0.1~1000mAs
- Input Power: 380V, three phases
- Support ABS function
- Support AEC function

X-Ray Tube

Specifications

- Brand: Siemens
- Model: SV150/40/80C-100_1
- Focus: 0.6mm/1.0mm
- X-ray machine
- Maximum output voltage: 150kV
- Nominal pulse fluoroscopy tube voltage: 110kV
- Tube anode heat capacity: 600kHU
- Target angle: 12°

Dynamic Flat Panel Detector HCP Version

Specifications

- Brand: Angell
- Detector type: flat panel, scintillator: CsI
- Imaging area: 43x43cm²
- Static pixel matrix: 3072*3072 (9.437 million pixels)
- Maximum Continuous Spot-Film pixel matrix: 1536*1536 (2.359 million pixels)
- Frame rate of Pulsed fluoroscopy: up to 12f/s
- Frame rate of Continuous fluoroscopy: up to 30f/s
- Spot-Film Preparation Time: ≤0.8s
- Output Grayscale: 16bit
- Imaging time: ≤6S, preview time ≤3S
- Static spatial resolution: 3.5 lp/mm

Remote Control Table

Specifications

- Brand: DC-DT-4
- Model: Angell
- Table Size: 2100 x 880 x 715mm
- Table Lateral Movement: 250mm
- SID: 1.15m~1.8m
- Table Rotate Range: -15° - +90°
- Tube Column Rotate Range: -35° - +35°
- Tube Assembly Swiveling Range: -90° , 0° , +90°
- Tube Arm Lateral Movement: 1300mm
- Tube Arm Vertical Movement: 650mm

- High Voltage Cables of X-ray Machine: Length: 8m
- Load Bearing of Table: 200kg
- Rotation Angle of Foot Pedal: 360°
- One-click in Place Function: Standing Position or Lying Position

Remote-Control Desk

Specifications

- Brand: YKT-4
- Model: Angell
- Desk Size: 1700 x 800 x 721 mm
- One-click in Place Function: One-click to Standing Position or Lying Position
- Exposure Fluoroscopy Switch: Exposure and fluoroscopy can be prohibited to ensure the safety of doctors and prevent accidental exposure

Movable Console with Monitor (Optional)

Specifications

- Brand: YKT-4
- Model: Angell
- Height Range: 770mm — 1170mm
- One-click in Place Function: One-click to Standing Position or Lying Position
- Exposure Fluoroscopy Switch: Exposure and fluoroscopy can be prohibited to ensure the safety of doctors and prevent accidental exposure

Collimator

Specifications

- Model: XS-F
- Brand: Angell
- Power: 150W(24VAC)
- Position Feedback of Operation: Yes
- Touch Display: Yes
- Filter: $\leq 2.0\text{mmAl@120KV}$

Grid

Specifications

- Brand: JPI

- Model: JPI GRID—1000
- Ratio: 10:1
- Focal Distance: 120cm

DX Ray Diagnostic Workstation

Specifications

- Brand: Angell
- Model: TS—5
- CPU: INTEL
- RAM: 8G
- Hard Disk: 1T
- Pixel matrix: 1280*1024 pixels
- Patient management: Manual registration, Automatic Query from Worklist
- Image acquisition: Automatic Window Adjustment, Automatic Sending, Static Image Acquisition, Dynamic Image Acquisition, Video Save, Playback, Bone Stitching
- Image processing: Image Correction, Image Clipping, Automatic Image Segmentation, EAE Image Enhancement Processing, IEQ Image Processing
- Image viewer: Window/Level Adjustment, Image Flip, Image Rotation, Image Magnification, Restore
- Medical report: Automatic Loading of Patient Information, Expert Template
- Film printing: Supports Standard DICOM 3.0 Printer
- DICOM transmission: Complies with the standard DICOM 3.0 PACS and Workstation
- Fluoroscopy amplification: 3 step, 17/15, 17/12, 17/9.
- Anatomically programmed radiography (APR) : auto pre—set parameter link to the chosen body part, at least 1350 preprograms.
- Automatic collimator: After choosing the body part, collimator can adjust the suitable window size, changing SID could also bring automatic adjustment of collimator

TOMO Wide range linear scanning depth Imaging (Optional)

According to the later clinical needs, it can be upgraded to large—scale linear scanning depth imaging (provided with third—party certification)

- Minimum layer thickness: 1mm (provide software screenshots for proof)
- Minimum spacing: 1mm (provide software screenshots for proof)
- Maximum number of layers: 300 layers (provide software screenshots to prove)
- Maximum coverage of a single imaging: 1100 mm (provide software screenshots as proof)

Shenzhen Angell Technology Co., Ltd.

Room 408, Building A, Huahan Innovatiosn Park,
Langshan Road, Nanshan District,
Shenzhen 518057, China

Tel: (+86) 755-86966744

Fax: (+86) 755-26859389

Web: www.szangell.com